



Manje vremena, manje krvi!

Element[®]

Stvara vrhunsku učinkovitost
sa minimalnim uzorkom



TUZLA-FARM d.o.o.
Rudarska 71
75000 Tuzla
Bosna i Hercegovina


Telefoni:
+387 35 321 202
+387 35 321 203
+387 35 321 204

PJ Sarajevofarm
Hasana Merdžanovića bb
71320 Vogošća
Bosna i Hercegovina

Telefon:
+387 33 660 711

Tačnost

Dijagram razlika tačnost sistema
(kombinacija 3 lota)



Specifikacija

Mjerenje	Kapilarna puna krv
Uzorak krvi	0,3 μl
Vrijeme mjerenja	3 sekunde
Raspon mjerenja	10~600mg/dL (0.6~33.3mmol/L)
Radna temperatura	10~40 °C (50~104 °F)
Radna vlažnost	10~90%
Hematokrit	20~60%
Memorija	365 EA
Rok trajanja baterije	Više od 5.000 kontinuiranih mjerenja
Veličina	82 x 55 x 18mm (H x W x D)

Karakteristike



Jednostavno korištenje

